



AGRICONS PP/PE CONTAINERS KIT

Disposable Sterile Medical Devices
(Class 2A CE0476 marked)

The screw cap containers have a double internal security closure and a translucent high resistance polypropylene body. They are designed for solid and liquid biological samples.

The containers are available in the following volumes:

- 20mL
- 60mL
- 120mL
- 500mL
- 40mL
- 80mL
- 250mL
- 1000mL

Disposable Sterile Medical Device
for Biological Materials Class 2A CE0476

Features

Cod. ADM1XXX

It's a container manufactured in translucent polypropylene (PP) with a resistance structure and a polyethylene (PE) screw cap with a double internal security closure;
It's double/triple packaging in a sterile double layer barrier:
BoPET/PE Film + Tyvek® 1073B 100% HDPE; Universal label;
Sterilization indicator;
Removable Label with double sterilization indicator;
Sterilization by irradiation accelerated electrons Beta;
Medical Device Class 2A CE0476.

Description

Cod. ADM1XXX

Usable for sampling, transportation and storage of biological materials.
Before using check the integrity of the package.
Respecting the proper operating rules in the asepse, insert the tissues inside the device by applying the cap and rotating it until it is completely closed.
In accordance with the internal procedures, properly transport or store the tissues, avoiding any possible thermal or mechanical shocks.
Following the proper operating rules in asepse, remove the tissue from the inside of the device by turning the cap up to a full opening.

AGRICONS RICERCHE Medical Device Division